

CoaguChek Pro II PT Tests



110678

UOM: 2X24

● Available to Order

CoaguChek PT Test (2 x 24) for Pro II

CoaguChek PT Test Strips are simple and easy to use, requiring a minimal blood sample in order to obtain a result. The sample of blood can be applied to the CoaguChek PT test strips on either side. These test strips are for use with the CoaguChek Pro II Monitor.

Unit of Measure: 2x24

This is a special order item and is non returnable, for an estimated delivery time please call 01323 445001.



Save As PDF

Demo Version - Select.Pdf SDK

Overview

Product Reviews

Sensitivity to heparin

- PT/INR: Insensitive to unfractionated and fractionated heparin concentrations up to 3 IU/mL blood
- aPTT: Insensitive to low molecular weight heparins (LMWH) up to 0.5 IU/mL blood

Quality control

On each strip, through the same channel as the blood passes

Stability

Store at + 2° C to + 30° C. Test strips can be used until expiry date printed on the box and test strip vial.

Box of 2 x tubes of 24 prothrombin time (PT) test strips for use with the CoaguChek Pro II PT/INR Monitor.