





Invacare Universal High Sling



- Support type:**
- Whole body and head
- Typical use:**
- Transfers from supine
 - Seated to seated transfers
- User characteristics:**
- Limited head, trunk and hip control

Size	A	B	C	D	Product number / Material		
					Polyester	Net	Spacer
XS	1040	740	610	340	1483776	1485950	1485955
S	1080	840	730	380	1483777	1485951	1485956
M	1420	1000	830	400	1483778	1485952	1485957
L	1550	1100	920	510	1483779	1485953	1485958
XL	1680	1150	920	550	1485949	1485954	1485959


Invacare Universal Standard Sling



- Support type:**
- Whole body including shoulders
- Typical use:**
- Seated to seated transfers
 - Good all round sling for the elderly
- User characteristics:**
- Only suitable for those with good head control

Size	A	B	C	D	Product number / Material		
					Polyester	Net	Spacer
S	1000	920	510	410	1485882	1478411	1485363
M	1090	1050	560	480	1485883	1478412	1485364
L	1180	1090	580	560	1485884	1478413	1485365
XL	1240	1250	585	710	1485885	1478414	1485366


Invacare Easy-fit Sling



- Support type:**
- Whole body and head
- Typical use:**
- Transfers from supine
 - Seated to seated transfers
- User characteristics:**
- Limited head, trunk and hip control
 - Amputees
 - Severely disabled

Size	A	B	C	D	Product number / Material	
					Polyester	Spacer
XS	1100	710	595	170	1451101	1451106
S	1130	760	640	180	1451102	1451107
M	1200	800	675	190	1451103	1451108
L	1270	860	720	210	1451104	1451109
XL	1390	940	780	220	1451105	1451110


Invacare Universal Low Sling



- Support type:**
- Hips and lower back
- Typical use:**
- Seated to seated transfers
- User characteristics:**
- Some trunk control but limited head control

Size	A	B	C	D	Product number / Material	
					Polyester	Spacer
S	850	800	350	410	1491893	
M	920	940	375	480	1491894	
L	950	1075	405	575	1451067	

Invacare Dress Toileting Sling (with or without head support)



- Support type:**
- Whole body and head (with head support)
 - Hips and lower back (without head support)
- Typical use:**
- Seated to seated transfers
- User characteristics:**
- Good head, trunk and hip control (without head support)
 - Limited head, trunk and hip control (with head support)

Size	A	B	C	D	Product number / Material	
					Toileting polyester	Toileting* polyester
S	870	710	200	430	1478440	1478444
M	950	800	230	490	1478441	1478445
L	1020	900	230	600	1478442	1478446
XL	1100	1080	270	760	1478443	1484780

* Version with added back support not shown




Yes, you can.®



Stand Assist Slings

Invacare Stand Assist Sling



Support type:

- Lower back

Typical use:

- Seated to standing transfers

User characteristics:

- Good body and head control
- Some weight bearing ability required
- Ability to grip boom

				Product number / Material	
Size	A	B	C	Polyester	
S	260	840	950	1484685	
M	330	895	1170	1484686	
L	360	950	1270	1484687	
XL	360	1020	1810	1484688	

Invacare Comfort Sling and Comfort Sling with Toilet Opening

Support type:


- Whole body and head

Typical use:

- Transfers from supine

User characteristics:

- Limited body control
- Amputee (above the knee)
- Chronic pain



					Product number / Material			
					Comfort		Comfort Toileting	
Size	A	B	C	D	Net	Spacer	Poly-ester	Net
XS	1040	970	510	960	1485967	1485948	1451091	1451096
S	1150	1050	550	1020	1485968	1485554	1451092	1451097
M	1400	1170	610	1180	1485969	1485777	1451093	1451098
L	1450	1290	670	1280	1485970	1485778	1451094	1451099
XL	1460	1330	750	1470	1485970	1451090	1451095	1451100

This sling offers a slightly reclined transfer position.

Invacare Transfer Stand Assist Sling

Support type:


- Lower back and hips

Typical use:

- Seated to seated transfers

User characteristics:

- Good body and head control
- Ability to grip boom



				Product number / Material	
Size	A	B	C	Polyester	
S	860	840	240	1484689	
M	1020	895	260	1484690	
L	1180	950	330	1484691	
XL	1230	1020	350	1485117	

Invacare Compatibility Statement

Invacare use a "Loop and Coat Hanger Spreader Bar System" as do many other manufacturers, therefore Invacare approves the use of its range of slings on suitable patient transfer systems manufactured by other companies. However, we do not recommend our slings be used on equipment other than the "Loop and Coat Hanger Spreader Bar System". In particular, it should be noted that **Invacare slings** are not suitable or compatible with "Keyhole Spreader Bar" or "Tilting Frame Spreader Bar" designs.

Sling Types	Fabric Matrix		
	Polyester	Net	Spacer
Universal Low	●		
Universal Standard	●	●	●
Universal High	●	●	●
Universal High Plus	●		
Dress Toileting	●		
Dress Toileting High	●		
Easy-fit	●		●
Amputee	●		
Comfort		●	●
Comfort Toilet	●	●	
Stand Assist	●		
Transfer Stand Assist	●		

Invacare Limited
 Pencoed Technology Park - Pencoed - Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1 656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail: ordersuk@invacare.com
 www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County Dublin - Ireland
 Tel: +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
 www.invacare.ie

© 2011 Invacare International Sàrl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare reserves the right to alter product specifications without prior consultation.

Sling Range - UK - 12/2011 - 1542214



Yes, you can.®

Size guide

