

CoaguChek Pro II System



110679

UOM: EACH

● Available to Order

The CoaguChek Pro II Coagulation Meter features multiple enhancements to the previous generation of CoaguChek devices

Unit of Measure: Each

This is a special order item and is non returnable, for an estimated delivery time please call 01323 445001.



Demo Version - Select.Pdf SDK

Overview

Product Reviews

It brings the lab closer to the patient in the hospital environment. In addition, it continues to support Vitamin K Antagonist (VKA) therapy monitoring in the professional segments, such as anticoagulation (AC) clinics and doctor's offices.

Two assays can be determined with the new CoaguChek Pro II meter - Prothrombin Time (PT) (with extended claims), as well as Activated Partial Thromboplastin Time (aPTT).

The CoaguChek Pro II meter features

2-D barcode reader for patient, operator and controls identification for faster ID entry

Universal battery pack, charging via power adapter and/or Base Unit *(HBU)

Enhanced connectivity through wireless technology and/or the unique QR code feature

PT/INR results in approximately 1 minute

aPTT results in less than 5 minutes

Side and top dosing options with a drop of venous, capillary or arterial fresh whole blood (min. 8 µL)

5 year warranty

Clinically Vital

From the floors to the OR to the ER to the ICU, aPTT and PT results are in your hands in real time

Enhanced Connectivity

Critical information is available when and where it is needed most

Wireless connectivity and bar code reader for digital input and output

QR code feature enables data transmission to any data management system immediately, with no need to transcribe results for a streamlined workflow