


TreadNought Floor Sensor Pad with Mono 1/4" Jack Plug



103195

UOM: Each

● In Stock

The TreadNought Nurse Call Mat works with the majority of installed Nurse Call Systems. The mat senses when a patient has stepped out of bed and triggers an alarm at the Nurse Call Station.

- Works via existing nurse call system
- Saves time - Reduces time-consuming room checks
- No false alarms, only activated when stepped on
- Increase patient safety
- Made from high quality double sided non slip grey rubber
- Strong welded seams

Unit of Measure: Each

Demo Version - Select.Pdf SDK



 Save As PDF

Overview

Product Reviews

Technical Specifications

- Dimensions
 - Width 60cm
 - Length 90cm
 - Depth 0.7cm
- Weight - 1.5kgs
- Pressure required to trigger the sensor - 2kgs

Spare leads are 110699 - Mono 1/4" Jack