

MobiCare

Anti-decubitus air cushion system designed for mobility
Compatible with any Carilex's power unit
Favorable for patients in sitting position

"MobiCare" is the mobilized air alternating pressure relief cushion designed to provide the same therapeutic function on the go without compromise. Besides MobiCare's standard equipped power unit, it is also connectible with any other Carilex power units to perform therapy for pressure ulcer prevention and treatment while patients need the therapy from bed to chair. MobiCare is a great accessory to go with any Carilex anti-decubitus air mattress system to provide patients a full range of choices under any environments.




carilex



carilex[®]

Power Unit Control Panel




Power Indicator
Power switch lights
when power is on

Comfort Knob



Foam base: Integrated 2"
foam base provides
additional protection in case
of power outage



Skid proof pad: Avoid seating
cushion from possible sliding



Quick release coupling: Easy
connection between power unit
and seating cushion

Technical Data

Individual Model No.

Power Unit: SR303
Cushion Set: MSR213

Electrical Input

AC100 – 120V/ 50 & 60 Hz,
or AC220 – 240V/ 50 Hz

Agency Approval

EN60601-1/ EN60601-1-2/ EN61000-3-2/
EN61000-3-3/ UL2601-1

Dimension

Power Unit: 26 x 22 x 8 cm
Cushion Set: 45 x 45 x 9 cm

Weight

Power Unit: approx. 2.5 kg
Cushion Set: approx. 1.0 kg

Power Consumption

Max. 10W

Material

Power Unit: ABS fire retardant UL
Cushion Cell: 100% nylon with TPU
lamination
Cushion Cover: 100% nylon with PU
backing

Operating & Transportation & Storage Temperature/ Humidity

5°C - 45°C / 15% - 60%

Maximum Weight Capacity

100 kg

This power unit is not AP/ APG
protected

carilex[®]

Turmbergstr. 9A, D-76227 Karlsruhe, Germany • Tel: +49(0)721 16 19 835 • Fax: +49(0)721 16 19 836
Email: info@carilex.de • Website: <http://www.carilex.de>